



Envi Elite Smart Chart



Envi Elite Smart Chart

The Envi Elite Acuity Chart comes fully equipped with all the necessary components for a standard acuity chart. The comprehensive kit comprises a Bluetooth remote, an HDMI device with software, a wall mount, and a high-definition monitor. The Envi Elite Acuity Chart is versatile, designed for use in both direct and mirrored exam rooms. Calibration is swift and user-friendly, providing expedited access to all charts and features.

Features:

- Contrast
- Cross cylinder
- Tumbling E
- Allen symbols
- Numbers
- ETDRS
- HOTV
- Landolt C
- Astigmatism Wheel
- Duochrome screen (Red/Green)
- Sloan or Snellen Optotypes (20/400 - 20/10)
- Mirrored lane compatibility
- Randomized optotype presentation
- Import videos and images
- Ishihara screening (Color test)
- Single line, single letter, or groups of lines and letters
- 1, 3, or 5 Optotypes per line
- Worth 4-Dot
- Metric/Customary (6/6 or 20/20)

