



envieye.com  
727-483-9140



# Envi HSL-20



## Envi HSL-20

The Envi HSL-20 is an advanced slit lamp microscope designed for high performance, accurate diagnostics, reliability, and patient comfort. High-quality optics result in excellent image quality and resolution. The HSL-20 features light tower LED illumination and two-step magnification.

### Technical Specifications:

Type	13 degree binocular tube with converging optics
Magnification Style	2 steps by objective lens rotation
Magnifications	10x, 16x
Eyepiece	10x
PD Adjustment	52mm-78mm
Diopter Adjustment	-8D - +5D
Slit Width	Continuous from 0mm to 14mm
Slit Height	Continuous from 1mm to 14mm
Aperture Diameter	14mm, 10mm, 5mm, 3mm, 1mm, 0.2mm
Slit Angle	0-180 degrees
Filters	Heat-absorbing, ND Filter, Red Free, Cobalt Blue
Slit Inclination	5, 10, 15, and 20 degree
Longitudinal Movement	110mm
Lateral Movement	110mm
Fine Base Movement	15mm
Vertical Movement	30mm
Fixation Bulb	LED

