

Huvitz

Distributed by:



AUTO-PHOROPTERS

HDR-9000

A cutting-edge phoropter designed for personalized and efficient refraction testing. Utilizing the innovative 21-point exam method, it simplifies the refraction process, making it accessible to all users. With its rapid astigmatism measurement, customizable settings, and wireless connectivity, the DR-9E delivers a superior examination experience.



FEATURES:



- **Simplified Refraction with 21-Point Precision:**
 - Simplifies refraction with the 21-point exam method, eliminating the need for complex knowledge or experience.
 - Displays clear results on the display for easy interpretation by both examiners and patients.
- **Rapid and Precise Astigmatism Assessment:**
 - Employs dual cross-cylinder and Jackson cross-cylinder lenses for highly accurate astigmatism assessments.
 - Improved lens movement speed minimizes accommodation interference, ensuring precise measurements.
- **Tailored Patient Examinations:**
 - Offers a tiltable body (0° to 45°) for comfortable near-vision examinations, simulating natural reading posture.
 - Provides monocular height adjustment (+/- 3mm) for personalized examinations in patients with varying eye heights.



CONTACT US

envieye.com | 833-627-0546
sales@nextvision2020.com