

envieye.com
727-483-9140



Envi LMA-20P



Envi LMA-20P

- Comes standard with Thermal Printer
- 5.6" full color touch screen that can be changed to black/white
- Wide range of measurement from -40D to +25D
- UV, Blue ray, and visible light transmittance measurement available
- Hartman Sensor enables measurement of a wide variety of lenses
- Auto calibration to calibrate fine errors during measurement
- Come standard with RS-232C and USB port for data transfer

Technical Specifications:

Sphere	-40.00D - +25.00D
Cylinder(+, ±, -)	0 - +10.00D
Axis	0° - 180° (1° Steps)
ADD	0 - 10.00D
Prism(X-Y/P-B)	0-10▲

