

envieye.com  
727-483-9140



envi

essential eyecare instruments



Envi S1-TC



## The Envi Family of Ophthalmic Exam Chairs

Each Envi chair was custom designed and built to serve the needs of the eye care provider based on the industries most used and requested features. The Envi line of chairs brings together unmatched **quality** and **affordability** in a contemporary package.

### Envi S1-TC | Tilt Chair

- Complete power chair (up/down)
- Chair up/down operated from side with soft-touch buttons
- Chair tilt control located on handle on back of the chair
- Adjustable headrest with easy lock and unlock
- Footswitch with up/down controls
- Chair rotates 180 degrees with foot lock on the backside of chair
- Folding footrest with spring assist
- High quality and durable upholstery
- Customizable upholstery colors available
- Chair base comes standard in white to match Envi S1-S Instrument stand



### Technical Specifications:

- Chair base: 21" x 23.75"
- Height in lowest position: 23.6"
- Height in highest position: 31.5"
- Length when fully reclined: 65"
- Lifting capacity: 500 lbs