

envieye.com
727-483-9140



Envi Switch



Envi Switch

The Envi Switch is a versatile stand that was custom designed and built to fit any practice size and tailors to the needs and wishes of any physician. The Switch is capable of the smallest footprint in the industry when paired with any Envi chair. It is also capable of serving as a stand-alone unit when fit into its custom baseplate. Switch it up with multiple configurations, unique add-ons, and practical storage.

Features:

- 3 charging wells and instrument tray
- Soft touch button controls for power of stand, slit lamp, overhead lamp, binding post, and two auxiliary inputs
- Control up/down of chair from slit lamp arm
- One hand operation of slit lamp arm
- Refractor arm with additional 3rd arm available
- LED Overhead lamp

Technical Specifications:

- Base size: 23.6" x 20.8"
- Total height of unit: 84.6"
- Slit lamp arm vertical range: 12"
- Slit lamp arm weight capacity: 80 lbs
- Phoropter vertical range: 10"

Custom baseplate to operate as traditional stand with functional storage spaces for lenses, drops, etc.





36" or 40" Footprint
← Right or Left handed setup →



envieye.com
727-483-9140



Multidirectional/dimmable LED overhead lamp with 360° rotation with lamp switch on slit lamp arm.

Technical Specifications:

- Base size: 12" x 24"
- Base size (with baseplate): 21.5" x 24"
- Total height of unit: 76"
- Slit lamp arm weight capacity: 80 lbs
- Phoropter vertical range: 10"



Paired with Envi S1-MC chair with secure anchoring rod for an edge-to-edge footprint of 36"