



# Envi Trial Lens Frame



# Envi Trial Lens Frame

A must-have for every eye care professional's toolkit. With its smooth controls and durable construction, it simplifies lens fitting, ensuring patient comfort and precise adjustments. Its familiar design and intuitive operation make it effortless to use.

Benefits include: comfortable fitting, lightweight frame, smooth and precise lens adjustment controls, long-lasting durability, and user-friendly operation.



## Features:

- Adjustable temple angle for personalized fit.
- Temple length adjustable from 4.125" to 5.5" (98 to 135mm) for enhanced comfort.
- Nose bridge height can be adjusted from 0 to 14mm.
- Nose angle can be adjusted from 0° to 30° for optimal positioning.
- PD (Pupillary Distance) adjustable from 24 to 40mm on each side (48mm to 80mm total for both eyes).
- Axial calibration on each lens:
  - Left lens: adjustable from 120° to 0° to 135°.
  - Right lens: adjustable from 45° to 180° to 60°.
- Lens rotation capability of full 360° for versatile usage.
- Accommodates three 38mm lenses in front and one behind for customizable vision correction.
- Clear, precisely printed scales for accurate adjustments.
- Rubber nose bridge for enhanced comfort during extended wear.

