

envieye.com
727-483-9140



Envi S1-AR



The Envi Family of Ophthalmic Exam Chairs

Each Envi chair was custom designed and built to serve the needs of the eye care provider based on the industries most used and requested features. The Envi line of chairs brings together unmatched **quality** and **affordability** in a contemporary package.

Envi S1-AR | Full Power Chair

- Folding footrest with spring assist
- Footswitch with up/down controls
- High quality and durable upholstery
- Customizable upholstery colors available
- Adjustable headrest with easy lock and unlock
- Armrests go up and down with chair when reclining
- Complete automatic power chair (up/down, recline back/forward, and auto)
- Chair rotates 180 degrees with foot lock on the backside of chair
- Chair base comes standard in white to match Envi S1-S Instrument stand
- Easily accessible soft-touch buttons on both sides of the chair



Technical Specifications:

- Chair base: 21" x 23.75"
- Height in lowest position: 23.6"
- Height in highest position: 31.5"
- Length when fully reclined: 65"
- Lifting capacity: 500 lbs