

envieye.com
727-483-9140



Envi Trial Lens Set




Envi Trial Lens Set

Discover the quality of the Envi Trial Lens Set, complete with Full-view lenses featuring silver polished metal rims for durability and clarity. Crafted with precision, this comprehensive set offers 266 lenses, including spherical, cylindrical, prismatic, and accessory variations. Built for reliability, it ensures consistent performance over time. Presented in an elegant carrying case, it allows for convenient storage and transport, serving as an essential tool for accurate vision testing in various clinical settings.

Features:

- Full-view construction for clear visibility and ease of examination.
- Versatile lens set covering a wide range of prescription needs, with a Sphere Range of (+/-) 0.25-20.0D and Cylinder of (+/-) 0.25-6.0D.
- Convenient organization within the carrying case for quick access to specific lenses.
- Compatibility with most trial frames for seamless integration into clinical practice.
- Enhanced patient comfort during vision testing with smooth, polished metal rims.



Trial Frame not included.