

# MARCO ENCORE AUTO CHAIR



## MARCO ENCORE AUTOMATIC

The Marco Encore Automatic Chair features an exclusive hydraulic system that controls both the elevation and recline motion. Easy-to-use backlit fingertip switches control all functions of the chair, as well as programmable recline positions, auto-recline/return, and safety functions. A corded foot pedal and single-handed operation headrest are also standard. Optional surgical accessories are available for minor surgical procedures.



EASY, "DIAL-IN"  
PROGRAMMING CONTROL  
WITH SAFETY SWITCH TO  
DISABLE CHAIR OPERATION

## ENCORE AUTOMATIC EXAMINATION CHAIR



Base: 22.25" x 23.85"

Lift Capacity: Approx. 375 lbs

Ground to Seat Height: 21" low – 32" high

Rotation: 270° degrees

Recline Mechanism: Auto Lay Back

Input Voltage: 120VAC 60HZ

Fuses: 8A 250V

Available Colors: Black, Charcoal, Navy, Silver, Silver & Black two tone, Silver & Navy two tone