NEXT
VISION INSTRUMENTS



iLight
LLLT Pro

Optometric Aesthetics

iLight LLLT Pro

Low-Level Light Therapy



LLLT, also known as photobiomodulation, utilizes low-level light to activate cellular processes and facilitate healing.

Benefit from direct optometrist-to-optometrist interaction, gaining practical, in-office instruction, as well as continual assistance and advice from experienced ODs.

Features:

- 2 Heads Included: All Red Head & Multi-Colored LED Head
- NO consumables. Heads last 35,000 hour before needing to be replaced
- Industrial wheels for easy movement to and from treatment rooms
- Extendable arm to accommodate patients sitting up, reclining, or lying flat for treatments
- Treats the full face and neck

$29,900 | 2-Year Warranty

All RED LED Head (633 nm):

- Irradiance up to 155mW/cm2
- Can heat facial tissue up to 43° celsius

Known uses in eyecare: OSD, MGD, Chalazion, Hordeolum, Blepharitis, Sjögens's disease, general inflammation, post procedure (bump removal, blepharoplasty, LASIK)

Red 633nm:

- Used for skin healing, improved texture, dry skin, post IPL/microneedling/chemical peels

Known uses in eyecare: (same as red light above but at a lesser energy level.

Blue 417nm:

- Used for acne, blackheads, redness/irritation, pigmentation, melasma, etc.

Known uses in eyecare: blepharitis, demodex, rosacea

Yellow 590nm:

- Used for fine lines, wrinkles, skin texture, tightening skin, etc

Known uses in eyecare: OSD, post procedure (bump removal, blepharoplasty, LASIK)

Infrared 850nm:

- Used to reduce swelling, post procedure (chemical peels, microneedling, injectables, IPL)

Known uses in eyecare: post procedure (bump removal, blepharoplasty, LASIK, cataract surgery)

Green 532nm

- Used to reduce inflammation, treat migraines and chronic pain conditions

NEXT VISION INSTRUMENTS

Doctor to Doctor

The main objective is to maintain a **direct Optometrist-to-Optometrist** experience, ensuring a seamless journey from exploring our devices to receiving hands-on, in-office training, all supported by *ongoing guidance from ODs*. Your education and mentorship will be overseen by a seasoned Optometrist.

Our efficient business model eliminates unnecessary layers to offer affordable medical devices that deliver impressive results. These devices not only benefit your patients but also **significantly boost your practice revenue,** ensuring a substantial return on investment for every Optometrist.

OA
Optometric Aesthetics